


Skin Roller Principle

Derma Rolling System roller has 192 (8 lines of needles) or 72 (3 lines of needles) micro needles that puncture the surface of epidermis, which blocks the active ingredients from entering into the skin. The skin roller will allow the active ingredients in the skin care products to absorb deep into the skin.



After skin needling, the penetration of skin care products active ingredients increases significantly when compared to the regular topical application. Due to increased skin care products absorption the skin gets a better dosage of highly recommended vitamin A and vitamin C contained in the skin care products, which helps the skin to release more collagen in a natural way.

Skin Rolling Benefits

Skin needling with the Derma Rolling System is an effective alternative for many skin care therapies. Main advantages of needling with the amazing skin roller are following:

- No skin damage.
- No fear of hyper pigmentation.
- Suitable for any part of the body.
- Suitable for all skin types.
- Suitable for thin and previously laser treated skin.
- Cost-effective compared to other skin therapies.
- Easy to learn and use.

Skin Needling Effects

Skin needling with our Make-up-Equipments skin roller induces skin's natural collagen building, that will help to restore skin tightness in the early stages of ageing, reduce pigmentation marks, thicken thin skin, smooth fine lines, soften different types of scars and reduce dilated blood vessels. The increased absorption of skin care products active ingredients will dramatically increase their effectiveness.

Skin needling with the Derma Rolling System has been used with success to:

- **Reduce scars** so after a few treatments they will flatten and become less obvious,
- **Reduce wrinkles**- lower-eyelid wrinkles, upper-lip lines and other facial wrinkles,
- **Reduce cellulite** for having beautiful thighs without any dimpling appearance,
- **Reduce stretch marks** so they almost fade.

Effectiveness of skin needling varies from one individual to the other. Most will see up to 50% of improvement, while others may get only up to 20% improvement. The most important is that using the skin roller is a safe way to improve your skin's appearance by enchasing its natural processes.

Instructions before use

We kindly request you to read the following instructions attentively. It contains important using information of Micro-needle Nurse System(MNS).

Verifying before use:

The roller has been sterilized and it should be delivered in a sealed clean bag and tube-container, so packing checking upon receipt are needed and do not use when there are packing damage. No return and no exchange after unboxing the roller.

What is the Micro-needle Nurse System(MNS)?

Traditional methods of nutrient can only infiltrate into the horny layer of epidermal, the effect is not suitable (only reach about 0.3% of effectiveness). To solve this problem, after numerous researches, scientists have developed a micro-needle therapy. This treatment can effectively increase the penetration capacity of nutrients through the epidermis and dermis of skin cells, is a new and miracle skin care technology for its maximal delivery method. The micro-needle therapy can stimulate dermal wound healing ability, thereby stimulate the skin and extend proliferation, make the skin epidermis proliferate about 8%, its effects can be comparable to the laser and scaled cosmetic surgery. With minimum side effects, safe and reliable, is the best choice for anti-aging treatment.

Theory of Micro-needle Therapy

Micro-needle therapy uses a kind of indented with 192 slender steel needles, accompanying with functional nutrient liquid, regularly rolling in the places which need to care, it can puncture more than 250,000 micro-needs on the epidermal within five minutes, form a efficient nutrient delivery system in the epidermis and subcutaneous tissue, with which the active ingredient of the nutrient can entry into skin effectively. Thus make skin maintenance, skin care, skin treatment and skin beauty to be perfect.

Micro-needle therapy has the advantages of safe, reliable and multi-function. It can be used to beauty, Body slimming, Hair regeneration, etc. It comes to be the preferred method for skin anti-aging.

Modes of Micro-needles Nursing Instrument:

| Models | 312-01 | 312-02 | 312-03 | 312-04 | 312-05 | 312-06 | 312-07 |
|---------------------------------|--------|--------|--------|--------|--------|--------|--------|
| Micro-needle depth | 0.25mm | 0.5mm | 0.75mm | 1.0mm | 1.5mm | 2.0mm | 2.5mm |
| Roller length/diameter | ● | ● | ● | ● | ● | ● | ● |
| Complete set (equipment) | ● | ● | ● | ● | ● | ● | ● |
| Micro-needle Tube | ● | ● | ● | ● | ● | ● | ● |
| Roller Cover cap | ● | ● | ● | ● | ● | ● | ● |
| Roller accessories (roller bag) | ● | ● | ● | ● | ● | ● | ● |
| Roller Bag | ● | ● | ● | ● | ● | ● | ● |
| Roller Cleaning | ● | ● | ● | ● | ● | ● | ● |
| Preparation items | ● | ● | ● | ● | ● | ● | ● |
| Roller product | ● | ● | ● | ● | ● | ● | ● |

Operation of MNS

- 1 Choose suitable model according to need.
- 2 Use medical alcohol to disinfect position need used and the equipment.
- 3 Do not suitable functional nutrition on skin need to be cared.
- 4 Roll back and forth 4-5 times orderly and keep equally in speed and pressure to create more micro-veasels.
- 5 Take your time and concentration.
- 6 Move the roller four times at least in four directions (where need cared), horizontally, vertically, and diagonally right and left. (refer to Figure)
- 7 Do not nutrition or mask products on skin where be cared and pay attention to defending skin inflammation.

Microneedle Moving Directions Figure: